

# AIR-CRUTCH™



## Walk on Air™


Introducing the smarter alternative to painful traditional crutches. The Air-Crutch combines clinical precision with user-first comfort, featuring Airgonomic™ patented underarm cushions with innovative air bladders designed to disperse pressure where users feel it most. With adjustable sizing to fit users from 4'7" to 6'7" and support up to 300 lbs, Air-Crutches are built to help a wide range of users move with less pain and experience a more comfortable recovery. For you, that means stronger customer loyalty, premium showroom appeal, and simplified inventory with one solution that fits users from youth to tall adult — all while protecting margins with MAP pricing.



### Advanced engineering where standard crutches hurt most

#### All-Day Comfort

With 57.9% thicker padding and 21% more surface area than standard crutch pads, the patented Airgonomic™ underarm cushions disperse pressure with layered air bladder and PU foam — for unmatched comfort.




#### Less Pressure, More Support

Naturally angled soft, gel handles reduce hand and wrist strain, delivering 27.6% less pressure than standard crutches.

#### Cushioned Confidence in Every Step

The patented Flex Grip Air-Tip softens every step with innovative air bladder technology, while the water-dispersing treads are designed to improve stability and helps users maneuver through wet or uneven surfaces.






#### Item #

#### Description

#### UPC

#### Product Specs

MD10300BK	Air-Crutch, Onyx	 822383030968
MD10300BL	Air-Crutch, Sapphire	 822383030975
MD10300GR	Air-Crutch, Slate	 822383030982

Height Range	Adjustable for users 4'7" - 6'7"
Product Dimensions	2.25" (W) x 7.25" (L) x 39.75" - 62.75" (H)
Product Weight	4.875 lb
Weight Capacity	300 lb
Warranty	Limited Lifetime
Retail Box Dimensions	24.5" (L) x 7.75" (W) x 5.75" (H)

