

## MD by Drive Air-Go Walker

HCPCS:




FIGURE A




FIGURE B



FIGURE C

Differentiate your offerings with the Air-Go Walker--delivering comfort, ease, and freedom beyond traditional walkers. With oversized, cushioned Airgonomic handles, water-dispersing treaded wheels, and a padded mesh seat, it gives users both freedom of movement and confidence in every step.

- **Airgonomic Comfort:** Oversized hand grips are easy to grasp and feature integrated air bladders that cushion the hands.
- **Convenient Seating:** Integrated, open-and-go mesh seat allows for rest anytime.
- **Natural Stride:** Curved  $\frac{3}{4}$ -length seat design allows users to walk comfortably without interference.
- **Easy, Fast Fold:** The simple scissor fold design allows for compact storage and easy transportation.
- **Improved Traction & Stability:** Large, water-dispersing treaded front wheels grip the earth help maneuver through diverse terrain.
- **Fully Assembled Out of the Box:** Ready to use immediately, requiring no tools or setup for added convenience.

| PRODUCTS              | Description                                |
|-----------------------|--------------------------------------------|
| MD30100BK             | Air-Go Walker, Onyx                        |
| MD30100BL             | Air-Go Walker, Sapphire                    |
| SPECIFICATIONS        | Description                                |
| Warranty              | Limited Lifetime                           |
| Product Weight        | 8 lb                                       |
| Material              | Anodized Aluminum                          |
| Seat Dimensions       | 21.5" (W) x 8.8" (D)                       |
| Handle Height         | 33.1" to 36"                               |
| Weight Capacity       | 300 lb                                     |
| Height                | Offering adjustments for users 4'9" - 6'7" |
| Seat To Floor Height  | 20 - 23                                    |
| Dimensions Assembled  | 25.6 (W) x 25.2 - 26.9 (L) x 33.1 - 36 (H) |
| Folded Dimensions     | 7.5" (W) x 26.2" (L) x 20.5" (H)           |
| Retail Box Dimensions | 21.5 (L) x 7.1 (W) x 28.5 (H)              |
| Wheel Size            | 6 (D) x 1 (W)                              |