



Wound Care

Essity has a full line of wound care products to support your athletes on and off the field.

Portfolio overview

STRETCHABLE LIGHT-
WEIGHT FABRIC DRESSING



Leukoplast® elastic

TRANSPARENT DRESSING
WITH ABSORBENT PAD



Leukomed® T Plus

WIDE AREA FIXATION



Leukoplast® Hypafix®

ELASTIC GAUZE BANDAGE



Elastomull®

HYDROPHOBIC BACTERIA
BINDING POST OP DRESSING



Leukomed® Sorbact®

CONFORMING, ADHESIVE
FIXATION BANDAGE



Cover-Roll® Stretch

WOUND & EYE
IRRIGATION



SplashCap®



EyeCap®

Leukoplast® elastic

Fabric Adhesive Bandages

- Absorptive, non-adherent pad
- Adhesive surrounds the island pad to seal the affected area
- Contours to different areas on the body
- Variety of shapes and sizes



formerly known as Coverlet®

REF/BNR	Size	Description	UOM
		Strips	
7645617	¾ in x 3 in		100/box
7645606	1 in x 3 in		100/box
7645604	1 in x 3 in		300/box
		Spots	
7645616	7/8 in in round		100/box
7645613	1¼ in in oval		100/box
		Patches	
7645612	1½ in x 2 in		100/box
7645610	2 in x 3 in		50/box
7645605	4 in x 2¾ in		50/box
		Knuckles	
7645607	1½ in x 3 in		100/box
		Toe Shield (non-sterile)	
7645618	2½ in x 1⅜ in		100/box
		Digits	
7645608	Sm. finger tip		100/box
7645609	Lg. finger tip		50/box
		4-Wing	
7645611	3 in x 3 in		50/box
7646600	—	Eye Occluder REG (non-sterile)	30/box
7646700	—	Eye Occluder JR (non-sterile)	30/box
7253702	—	Cabinet (empty)	1 unit
7253703	—	Cabinet (filled)	1 unit

Leukomed® T Plus

Transparent, Total-Island Dressing with Absorbent Non-Adherent Pad

- Waterproof and barrier to bacteria
- Easy to apply with the unique red stripe system
- Should clinical conditions allow, the dressing may remain in place for up to 7 days




REF/BNR	Dressing Size	Wound Pad Size	UOM
7238200	2 in x 3 in	1 1/8 in x 1.5 in	50/box
7238201	3 1/8 in x 4 in	1.5 in x 2 3/8 in	50/box
7238202	3 1/8 in x 6 in	1.5 in x 4 3/8 in	50/box
7238214	4 in x 7 7/8 in	2 in x 6 3/8 in	50/box
7238203	4 in x 10 in	2 in x 8 1/8 in	50/box
7238204	4 in x 11 7/8 in	2 in x 9 3/4 in	50/box
7238205	4 in x 13 3/4 in	2 in x 11 3/4 in	50/box

4.3 WIDE AREA FIXATION

Leukoplast® Hypafix®

Wide Area Fixation

- Secure, rapid fixation that is excellent on contoured parts of the body, can be used as **Turf-Tape**
- Entire wound dressing is protected against contamination
- Convenient cut-to-size roll eliminates unnecessary waste



REF/BNR	Size	UOM
7193921	2 in x 2 yds (5 cm x 2 m)	36/box
7193919	4 in x 2 yds (10 cm x 2 m)	12/box
7193920	6 in x 2 yds (15 cm x 2 m)	12/box
7144301	2 in x 11 yds (5 cm x 10 m)	24/box
7144302	4 in x 11 yds (10 cm x 10 m)	24/box
7144303	6 in x 11 yds (15 cm x 10 m)	24/box

Elastomull®

Elastic Gauze Bandage

- High density fabric for a soft, rich texture
- Performance and stretchability to help control bleeding and edema
- Precision fit in areas with constant movement
- Helps prevent slippage once bandage is in place



REF/BNR	Size	UOM / Case
Non-sterile		
0209400	1 in x 4.4 yds	20 rolls
0209500	2 in x 4.4 yds	20 rolls
0209600	3 in x 4.4 yds	20 rolls

4.5 HYDROPHOBIC BACTERIA BINDING POST OP DRESSING

Leukomed® Sorbact®

Post-Operative Dressing with Sorbact® Technology

- Transparent and waterproof film with its wound contact layer coated with Sorbact® Technology
- Sorbact® Technology binds, inhibits, and removes hydrophobic bacteria utilizing a physical mode of action, without releasing active substances into the wound
- Reliable adhesion that is breathable and shower-proof
- Absorbent wound pad for low to moderate wound exudate
- Can be worn for up to 7 days depending on the clinical conditions of the wound
- For low to moderate exudate
- Sterile and not made with natural rubber latex



REF/BNR	Size	UOM
7619935	2 in x 3 in	50/box
7619936	3 in x 4 in	20/box
7619937	3 in x 6 in	20/box
7619938	4 in x 8 in	20/box
7619939	4 in x 10 in	20/box
7619940	4 in x 12 in	20/box
7619941	4 in x 14 in	20/box

Cover-Roll® Stretch

Conforming, Adhesive Fixation Bandage

- Stretchable non-woven polyester material adapts well to body contours, can be used as **Turf-Tape**
- Highly conformable to contour to all parts of the body requiring flexibility
- Secure fixation

REF/BNR	Size	UOM
4554700	2 in x 2 yds	36 rolls/case
4554800	4 in x 2 yds	12 rolls/case
4554900	6 in x 2 yds	12 rolls/case
4555000	8 in x 2 yds	12 rolls/case
4555100	12 in x 2 yds	12 rolls/case
4555200	2 in x 10 yds	12 rolls/case
4555300	4 in x 10 yds	12 rolls/case
4555400	6 in x 10 yds	12 rolls/case
4555500	8 in x 10 yds	12 rolls/case
4555600	12 in x 10 yds	12 rolls/case



4.7 WOUND & EYE IRRIGATION

SplashCap®

Irrigate and flush wounds rapidly

- Rapidly flushes contaminants from wounds
- Acute wounds
- Use with standard irrigation bottle not included

REF/BNR	Color	UOM
7446800	Clear	25/case
7446801	Clear	100/case



EyeCap®

Irrigate and flush eyes quickly

- Rapidly flushes contaminants from eyes
- Eye contamination
- Use with standard irrigation bottle not included

REF/BNR	Color	UOM
20025	Clear	25/case

